

Sci-70 ANTIGEN

AROTEC_Sci-70_Product_Info.pdf Version/Date: B/04.05.20

ATS01-02	Sci-70 antigen	0.20 mg
ATS01-05	Sci-70 antigen	0.50 mg
ATS01-10	Sci-70 antigen	1.0 mg

Description of the Product

Purified from bovine thymus. After coating onto ELISA plates the product will bind autoantibodies to Sci-70.

Purity: The Sci-70 autoantigen (75-82 kDa) is more than 90% pure, as assessed by SDS gel electrophoresis.

Concentration: 0.1-1.0 mg protein/ml.

Storage: The product is stabilised with 20% glycerol and 0.1% Mic-O-protect. Store at -20 °C or below (long term) or at +4°C (short term). Avoid repeated freezing and thawing. Mix thoroughly before use.

Clinical and Biochemical Data

Scleroderma (progressive systemic sclerosis) is a systemic autoimmune disease characterised by collagen deposition and connective tissue destruction of the skin, blood vessels and certain internal organs¹. The prevalence of scleroderma has been estimated to be 1 per 5000². Approximately 20-28% of scleroderma patients have autoantibodies to a nuclear protein referred to as Sci-70³. Anti-Sci-70 antibodies are found with a greater frequency in systemic sclerosis patients with proximal skin involvement than in systemic sclerosis patients with limited skin involvement⁴.

In 1986 the Sci-70 antigen was identified by Shero et al.⁵ to be the superhelical DNA-relaxing enzyme topoisomerase I. This report has since been confirmed by others^{6,7}. The molecular weight of enzymatically-active topoisomerase I determined by SDS-electrophoresis has been reported to range from 67-100 kDa⁸. Recent data suggest that the 100 kDa form is the predominant form *in vivo*, and that the lower molecular weight forms are due to proteolytic degradation⁹. Both the native and lower molecular weight forms have been shown to effectively bind anti-Sci-70 autoantibodies¹⁰.

The amino acid sequence of human Sci-70 antigen (DNA topoisomerase I) has recently been described after cloning and sequencing its cDNA¹¹. The deduced molecular weight of the mature protein is 90,649 kDa; lower than the SDS-electrophoretic molecular weight. This may be due to post-translational processing or to anomalous electrophoretic behaviour of the protein. The protein was found to be very basic, with a calculated isoelectric point of 10.05, and to possess its catalytic activity within a 67.7 kDa C-terminal fragment. Although the sequence of bovine Sci-70 has not been determined, the utility of the bovine antigen for the detection of autoantibodies to Sci-70 in patient sera has been demonstrated by several authors^{10, 12-14}.

Methodology

The following is an ELISA procedure which can be used to detect anti-Sci-70 autoantibodies in human serum using the ATS01 purified antigen:

1. Dilute the purified antigen to 0.5-1.0 µg/ml in PBS (10 mM potassium phosphate, pH 7.4, 0.15 M NaCl).
2. Coat ELISA plates with 100 µl of diluted antigen per well. Cover and incubate 24 hours at +4°C.
3. Empty the plates and remove excess liquid by tapping on a paper towel.

4. Block excess protein binding sites by adding 200 µl PBS containing 1% BSA per well. Cover and incubate at +4°C overnight.
5. Empty plates and apply 100 µl of serum samples diluted 1:100 in PBS / 1% BSA / 1% casein / 0.1% Tween[®] 20. Incubate at room temperature for 1 hour.
6. Empty plates and add 200 µl PBS / 0.1% Tween[®] 20 per well. Incubate 5 minutes then empty plates. Repeat this step twice.
7. Apply 100 µl anti-human IgG-enzyme conjugate (horseradish peroxidase or alkaline phosphatase) diluted in PBS / 1% BSA / 1% casein / 0.1% Tween[®] 20 per well and incubate for 1 hour.
8. Repeat step 6.
9. Add enzyme substrate and stop the reaction when appropriate.
10. Read absorbance in an ELISA spectrophotometer.

References

1. Medsger, T.A. Jr. (1985) in *Arthritis and Allied Conditions*, 10th Ed. (McCarty, D.J. ed.), pp 994-1036, Lea & Febiger, Philadelphia.
2. Steen, V.D. & Medsger, T.A. Jr. (1990) *Rheum. Dis. Clin. North Am.* **16**, 1
3. Douvas, A.S. et al. (1979) *J. Biol. Chem.* **254**, 10514
4. Steen, V.D. et al. (1988) *Arthritis Rheum.* **31**, 196
5. Shero, J.H. et al. (1986) *Science* **231**, 737
6. Guldner, H. et al. (1986) *Chromosoma* **94**, 132
7. Maul, G.G. et al. (1986) *Proc. Natl. Acad. Sci. USA* **83**, 5145
8. Schmitt, B. et al. (1984) *Eur. J. Biochem.* **144**, 127
9. Samuels, D.S. et al. (1994) *Mol. Biol. Reports* **19**, 99
10. Hildebrandt, S. et al. (1990) *Clin. Immunol. Immunopathol.* **57**, 399
11. D'Arpa, P. et al. (1988) *Proc. Natl. Acad. Sci. USA* **85**, 2543
12. Juarez, C. et al. (1988) *Arthritis Rheum.* **31**, 108
13. Tsay, G. et al. (1990) *J. Rheumatol.* **17**, 1314
14. Seelig, H.P. et al. (1993) *J. Immunol. Methods* **165**, 241

Mic-O-protect is from Roche Diagnostics GmbH (Mannheim, Germany).

Tween[®] 20 is a registered trademark of ICI Americas Inc.

NOTE: **No patented technology** has been used by AroTec during the preparation of this product.



AROTEC DIAGNOSTICS LIMITED