

Product information:

Product Name: Alternative names:	LC-1 antigen Liver Outosol type 1 antigen: formiminotransferase cyclodeaminase: ETCD
	Liver Cytosof type 1 antigen, forminiouralisierase cyclobeaninase, 1 rcb
Description:	Full length human sequence FTCD enzyme (isoform A) with hexa-Histidine tag
Source Material:	Recombinant human protein expressed in Spodoptera frugiperda Sf21 insect cells
Purification:	Purified using conventional chromatography procedures
Purity:	>90% pure (by SDS-PAGE)
Appearance:	Clear colourless solution

Antigen Overview:

LC-1 is a marker for autoimmune hepatitis (AIH) type II and antibodies against LC-1 are found to be detected in up to 50% of cases. Immunofluorescence of LC-1 is often not identified due to the brightness of LKM-1 obscuring accurate visualisation. Some occurrences of anti-LC-1 in conjunction with antinuclear antibodies (ANA) and smooth muscle antibodies (SMA) in patients with chronic hepatitis C and autoimmune hepatitis, is known. Some case studies have shown that anti-LC-1 is found as the only serological marker in patients with AIH type II.

The antigen recognised by LC-1 antibodies has been identified as formiminotransferase cyclodeaminase (FTCD). FTCD is a bifunctional enzyme which consists of 8 identical molecular subunits of 60 kDa each arranged in a circular tetramer of dimers. FTCD catalyses two independent but sequential reactions in the histidine degradation pathway in the mammalian liver. Each subunit consists of an N-terminal transferase active domain and a C-terminal deaminase active domain. The dimers are responsible for the separate activities of both the transferase and deaminase activity.

Ordering Information:

ATL04-02 - LC-1 antigen (recombinant) - 0.20 mg ATL04-10 - LC-1 antigen (recombinant) - 1.0 mg Custom pack sizes available on request

Storage Conditions & Handling:

Store at -65°C or below Avoid repeated freezing and thawing Storage buffer contains 20% Glycerol as cryoprotectant Mix before use and keep on ice

Applications:

After coating onto ELISA plates the product will bind auto antibodies to LC-1; recommended starting ELISA coating concentration 0.2 µg/well

Positive on Western Blot using sample with auto antibodies to LC-1 (e.g. patient serum), antibodies against FTCD or hexa-Histidine tag

Product Profile:



Fig. 1

Fig. 1. SDS PAGE of ATL04 LC-1 antigen; Coomassie stained; 2ug/lane Fig. 2. Western Blot of ATL04 LC-1 antigen with LC-1 patient serum (1:300 dilution). Colorimetric detection

Fig. 3. ELISA titration assay of ATL04 LC-1 antigen using LC-1 patient serum (1:100 dilution)

AROTEC_LC-1_Datasheet B/26.11.19

References:

MacKenzie RE, Aldridge M, Paguin J. (1980) The bifunctional enzyme formiminotransferase-cyclodeaminase is a tetramer of dimers. J. Biol. Chem. 255:200-211

Muratori L et al. (1995) Detection of anti-liver cytosol antibody type 1 (anti-LC1) by immunodiffusion, counterimmunoelectrophoresis and immunoblotting: comparison of different techniques. J. Immunol. Meth. 187: 259-63

Muratori L et al. (1998) Liver/kidney microsomal antibody type 1 and liver cytosol antibody type 1 concentrations in type 2 autoimmune hepatitis. Gut 42:721-6

Muratori L et al. (2001) Distinct epitopes on formiminotransferase cyclodeaminase induce autoimmune liver cytosol antibody type 1. Hepatology 34:494-501

Beland K et al. (2004) Anti-LC1 autoantibodies in patients with chronic hepatitis C virus infection. J. Autoimmun, 22:159-66

Lenzi M et al. (1995) Liver cystosolic 1 antigen-antibody system in type 2 autoimmune hepatitis and hepatitis C virus infection. Gut 36: 749-54

Murley LL, Mejia NR, MacKenzie RE (1993) The nucleotide sequence of porcine formiminotransferase cyclodeaminase. Expression and purification from Escherichia coli. J. Biol. Chem. 268: 22820-4

Kohls D et al. (2000) The crystal structure of the formiminotransferase domain of formiminotransferase cyclodeaminase: implications for substrate channeling in a bifunctional enzyme. Structure 8: 35-46

Lapierre P et al. (1999) Formiminotransferase cyclodeaminase is an organ-specific autoantigen recognized by sera of patients with autoimmune hepatitis. Gastroenterol. 116:643-9

Martini E et al. (1988) Antibody to liver cytosol (anti-LC1) in patients with autoimmune chronic hepatitis type 2. Hepatol. 8:1662-6

Renous R et al. (2004) Characterization of the antigenicity of the

formiminotransferase-cyclodeaminase in type 2 autoimmune hepatitis. Exp. Cell Res. 292: 332-41 Rigopoulou El et al. (2007) Autoimmune hepatitis-specific antibodies against soluble liver antigen and liver cytosol type 1 in patients with chronic viral hepatitis. J. Autoimmun. Dis. 4:2-7

Postal Address: PO Box 38586 Wellington Mail Centre Wellington 5045 New Zealand

Physical Address: 207 Gracefield Road Lower Hutt Wellington 5010 New Zealand

Contact Information: Phone: +64 4 569 0361 Fax: +64 4 569 0366 Email: info@arodia.com Web: www.arodia.com